

Shipping Specimen

For actual use please follow the directions in your kit.

STEP 3: Shipping your blood specimen card

ALLOW BLOOD SPOT TO AIR DRY FOR 30 MINUTES BEFORE BEGINNING THIS STEP.

Step 3.a. - Place the Blood Specimen Collection Card in Specimen Return Pouch after Drying



Fold the protective flap back to its original position on the Blood Specimen Collection Card. Place the Blood Specimen Collection Card in the Specimen Return Pouch.

Step 3.b. - Tightly Seal Specimen Return Pouch



Peel off the adhesive cover at the top of the Specimen Return Pouch. Fold down along the perforation. Press down hard along the back of the fold to seal it tightly.

Step 3.c. - Place Specimen Return Pouch in Cardboard Envelope



Place the sealed Specimen Return Pouch in the Cardboard Envelope. To seal the Cardboard Envelope, peel off the adhesive cover, fold down the flap and press down hard on the flap.

Step 3.d. - Place Cardboard Envelope in a FedEx® Shipping Envelope



Place the sealed Cardboard Envelope into the pre-paid Shipping Envelope. Seal the shipping envelope by peeling off the adhesive cover and folding the flap over.

Step 3.e. - Use FedEx® for Shipping

